

#### RX ONLY

# This product is intended only for use as described.

## What is Alzair™?

Alzair is a clinically proven nasal powder that works to strengthen your resistance to airborne allergens, such as pollen, dust, animal dander, and mold, which are inhaled through your nose. Alzair is a natural, plant-based extract that meets the highest purity and safety standards.

### How does Alzair™ work?

Studies have shown that by naturally enhancing nasal mucus, it is possible to filter out airborne allergens. Alzair works by forming a colorless, mucus-like gel lining in the nasal tract that acts as a barrier against pollen, dust, animal dander, and mold.

#### Indications for Use

Alzair Allergy Blocker is intended to treat hay fever and allergy sufferers by promoting alleviation of mild allergic symptoms (i.e. mild nasal irritation including itchy, runny, or congested nasal passages) triggered by the inhalation of various airborne allergens including indoor and outdoor environmental pollen, dust, animal dander, and mold.

Application of Alzair produces a protective gel-like barrier that evenly coats the nasal membranes and acts to block inhaled allergens within the nasal cavity.

Ingredients: Hydroxypropyl methylcellulose, mint scent

#### **Adverse Reactions**

Excessive mucoid or gel buildup

# Warnings

For nasal use only.

Do not use if tamper evident seal is broken or missing. Keep out of reach of children.

## Cleaning the Bottle

Should the bottle touch the nose or be contaminated by moisture, before replacing the cap, wipe the bottle thoroughly with a clean, dry, unused tissue. Only replace the cap when you are sure that no sign of dampness or contamination remains. Never use a bottle if you suspect it is contaminated.

Keep the bottle and especially nozzle dry at all times. Do not rinse the bottle with water or any other liquid as this may block it. Store at room temperature, away from moisture.

Once opened, use within 6 months

### Who can use Alzair™?

Alzair is suitable for adults, including pregnant or breastfeeding women, and for pediatrics\*.

#### How to use Alzair™

Before first use, test the pressure needed to administer an ideal dose. To do this, shake the bottle and remove the cap. Aim the bottle away from yourself and squeeze the bottle into the air.

For pediatric patients, squeeze gently for a 2-inch plume.


For patients 12 years and up, squeeze firmly for a 6-inch plume.

## If using a medicated nasal spray, use prior to Alzair

- 1. Gently blow your nose. Breathe out.
- 2. Shake the Alzair bottle.
- 3. Squeeze the bottle away from yourself to test the pressure needed to administer an ideal dose: 6-inches for patients 12 years and up, 2-inches for pediatric patients, .
- 4. Place a finger over one nostril to close it.
- 5. Place the Alzair bottle nozzle in the other nostril.
- 6. Quickly and firmly push the sides of the bottle together to deliver one "puff" of Alzair powder while gently inhaling.
- 7. Wait two seconds and then gently inhale to enable the Alzair powder to penetrate into the nasal passages. Repeat steps 3-6 on the other nostril.
- 8. Alzair may be used as often as required, but it is recommended three times a day (minimum). The key to getting the best out of Alzair is to make sure you maintain a constant layer of powder across the inner lining of your nose. Doing this creates the barrier between the aggravating allergens and sensitive membrane within the nasal tract.
- 9. Take Alzair before entering an environment containing airborne allergens
- e.g., going out to the garden or dusting.
- 10. Remember to reapply each time you blow your nose.



**RX ONLY** 











Manufactured For:

Rev. 01/18

594-599] of patients 8-18 years old with seasonal allergic rhinitis who were successfully treated

Based on the results of a clinical trial reported in Pediatric Allergy and Immunology [22 (2011);

Nasa ALL USA RX 08/17 50204 A

Hudson Scientific LLC Hudson, MI

For questions or information visit www.alzair.com or email info@alzair.com Call toll-free 888-292-8080